Attorney Docket No. V1025/20184 Customer No. 03000

PATENT

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE APPLICATION FOR LETTERS PATENT

5 APPLICANT

Robert M. Danen

POST OFFICE ADDRESS

1288 Avalon Drive

San Jose, California 95125

INVENTION

METHOD FOR CORRECTING VESSEL AND BACKGROUND LIGHT INTENSITIES USED IN BEER'S LAW FOR LIGHT SCATTERING IN TISSUE

10

15

ATTORNEYS : Caesar, Rivise, Bernstein,

Cohen & Pokotilow, Ltd.

12th Floor, Seven Penn Center

1635 Market Street

Philadelphia, PA 19103-2212

TO ALL WHOM IT MAY CONCERN:

20

Be it known that I, Robert M. Danen, a citizen of the United States, residing in the town of, County of Santa Clara, State of California, have made a certain new and useful invention in a METHOD FOR CORRECTING VESSEL AND BACKGROUND LIGHT INTENSITIES USED IN BEER'S LAW FOR LIGHT SCATTERING IN TISSUE, of which the following is a specification:

25

SPECIFICATION

FIELD OF THE INVENTION

The invention pertains to methods for intravital microscopy and more particularly to a method for correcting the vessel and background intensities used in Beer's Law for light scattering in tissue from an orthogonal polarized spectral imaging system.

BACKGROUND OF INVENTION

A recent innovation in the field of intravital microscopy is the introduction of a non-invasive instrument that produces high-contrast microvascular images which can be used in surgery, research and diagnostic applications. In particular, orthogonal polarization spectral (OPS) imaging systems (e.g., the device sold under the tradename CYTOSCAN® by Cytometrics LLC of Exton, Pennsylvania) utilize non-invasive instruments that produce an image by using scattered polarized light which creates a "virtual" light source within the observed tissue. The images generated by the OPS imaging systems allow doctors/clinicians the ability to visualize and measure real time images of the living being's microcirculation without the need for fluorescent dyes or transillumination.

5

10

15

20

However, as will be discussed in detail later, the photons from the light transmitted from the OPS imaging systems to the underlying tissue structure, travel in all directions. These photons scatter in all directions within the tissue, thereby distorting the relative intensities between blood vessel regions and background regions.

The following discussion relates to current references or activities associated with microvessel imaging.

"Modeling and simulation of illumination effects for evaluation of microvessels of the conjunctiva" (Wick, Loew, and Kurantsin-Mills 1996, *Modeling in Physiology*, H1229-H1239) involves modeling the reflection of source light within the conjunctiva.

To the best of Applicant's knowledge, researchers (Chance, et al., Yodh, et al.) at the University of Pennsylvania have developed techniques that measure oxygenation of hemoglobin. Jacques (Univ. of Oregon) has developed a polarization method to study the structure of epidermal tissue. However, the techniques referred to by Chance, et al.,

Yodh, et al. and Jacques measure optical properties of blood and tissues generally at near-infrared wavelengths. Most attempts to reconstruct both the absorption and scattering of tissue are on macroscopic scales (~0.5mm). At near-infrared wavelengths, light scattering dominates absorption; thus, elaborate inversion techniques are needed to determine the tissue optical parameters.

Thus, there remains a need for a method that solves for light scattering as a perturbative method in order to correct an absorption image when tissue vessels reside in the focal plane of an orthogonal polarization spectral (OPS) imaging system.

SUMMARY OF THE INVENTION

10

15

20

5

A method for correcting light intensities from blood vessels and background tissue beneath tissue surface in a living being that are exposed to polarized light from an orthogonal polarized spectral (OPS) imaging system. The method comprises the steps of: emitting polarized light at tissue comprising blood vessels and background tissue, wherein the blood vessels are located at a focal plane of the OPS imaging system and wherein a foreground region is formed between the focal plane and the tissue surface; collecting depolarized light that has impacted the blood vessels and that has experienced scattering within the foreground region and wherein the de-polarized light emerges from the tissue surface; generating a first image based on the collected de-polarized light; estimating the intensity of light scattered in the foreground region; and subtracting the intensity of light scattered in the foreground region from the first image to generate a corrected image based on focal plane light intensities.

DESCRIPTION OF THE DRAWINGS

Fig. 1 is a block diagram of an exemplary apparatus that uses orthogonal polarization spectral (OPS) imaging and depicts how polarized light, emitted from the OPS probe at target tissue, interacts with underlying tissue structure and is captured by the probe for processing;

Fig. 2 is a physical model of light propagation within tissue that originated from an orthogonal polarization spectral (OPS) imaging system;

Fig. 3 is a flow diagram of the method of the present invention; and

Fig. 4 is a flow diagram of the software of the present invention.

5

10

15

20

DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENTS

The invention is a method to correct the intensities of an absorption image for light scattered within foreground tissue using orthogonal polarization spectral (OPS) imaging. Specifically, the amount of light scattered between the focal plane of the probe and the tissue surface is estimated from the image. This scattered light contribution is then subtracted from the original image.

Before the method of the present invention is discussed, a summary of an exemplary device using OPS is described, such as the device marketed under the tradename CYTOSCAN® by Cytometrics LLC of Exton, Pennsylvania. As shown in Fig. 1, the device 1 comprises a probe 2 and a base unit 3. The base unit 3 comprises, among other things, a processor 4 (e.g., Pentium II processor), a PCI interface card 5, a display 6, a light source 7 (e.g., Welch-Allyn metal-halide arc lamp), an uninterruptible power supply (UPS) 8, and a keyboard 9. The probe 2 is electrically and optically coupled to the

base unit 3 through a harness (not shown).

5

10

15

20

The probe 2 comprises a lens assembly 10 that emits a beam of polarized light 11 filtered by filter 10A at 550 nm (the importance of which is discussed below) aimed at the target tissue 12. The majority (e.g., 90%) of the polarized light 11 is reflected from the target tissue surface 12A as glare 13A while the remaining portion of the polarized light penetrates deeply through the skin and impacts underlying tissue 12B, causing multiple scattering events 14. When this polarized light of 550 nm impacts various structures, it back-lights the structures, causes hemoglobin 12C to fluoresce (550 nm is the isobestic point for hemoglobin) and de-polarizes/scatters the light. This de-polarized light 13B is reflected back to and out of the target tissue surface 12A. As a result, the probe 2 receives both the glare 13A and the de-polarized light 13B which are passed through a lens 15 and an analyzer/orthogonal polarizer (APO) 16 inside the probe 2. The APO 16 (which is orthogonal to the first polarizer 10B) filters out the glare 13A and passes the de-polarized light 13B to a CCD video (or digital) camera 17 which converts the optical signal to an electrical (e.g., digital) signal for processing by the processor 4 and ultimately for display on the display 6 where an image (e.g., red blood cells, vessels, etc.,) is formed of the structures impacted by the incoming light and the scattered light.

When the tissue vessels (e.g., blood vessels, BV) are positioned in the focal plane F (Fig. 2) of the OPS probe 2, an incorrect estimation of the hemoglobin (Hb) occurs. In particular, light scatters in foreground tissue 32, the region of tissue between the tissue surface 12A and focal plane F of the OPS imaging system; a diffuse background 34 region is also shown in Fig. 2. Because the light field is diffuse in the focal plane F, photons travel in all directions, not just the direction defined by the collection optics (e.g., lens 15/AOP 16)

of the probe 2. These photons can scatter within the tissue 12B, distorting the relative intensities between vessel and background regions. This scattered light distorts the contrast of objects in the focal plane F. Because the concentration of Hb [Hb] is determined by measuring the quantity of light absorbed by vessels (e.g., blood vessels BV) in the focal plane F, light scattering in the foreground tissue 32 degrades measurements of [Hb].

Generally, by measuring the intensity of light absorbed by a blood vessel BV, [Hb] can be determined using Beer's law (A = 2-log₁₀%T, where A is absorbance and T is transmittance). However, Fig. 2 shows a diagram of light propagation within surface tissue for the OPS imaging system. As a result of the design of the OPS imaging system, all detected light travels from below the focal plane before reaching the tissue surface and the light field is diffuse within the focal plane F. The intensity of light measured by the OPS probe 2 at position x on the tissue surface ($I^{M}(x)$) is the sum of two components: light traveling directly from position x on the focal plane ($I_{f}(x)$) and light scattered from other positions on the focal plane F, ($I_{s}(x)$):

$$I^{M}(x) = I_{f}(x)e^{-\mu_{i}z_{f}} + I_{s}(x)$$
 , where

5

10

15

20

 $I^{M}(x)$ = light intensity measured at position x on the tissue surface 12A;

 $I_f(x)$ = light intensity that travels without scattering from position x on the focal plane (F) to position x on the tissue surface 12A;

 $I_s(x)$ = light intensity that scatters at 14 from other positions x_i ($x_i \neq x$) on the focal plane (F) to position x on the tissue surface 12A.

In the above equation, $\,\mu_t,\,$ is the transport coefficient of the scattering medium and z_f is

the distance between the focal plane F and tissue surface 12A.

5

10

15

20

As mentioned previously, use of measured intensities I^M as representative of the light distribution in the focal plane F where vessels reside, yields an incorrect estimation of [Hb]. Thus, a correct determination of [Hb] using Beer's Law requires that the light intensity distribution in the focal plane F be extracted.

Referring now in detail to the various figures of the drawing wherein like reference characters refer to like parts, there is shown at 20, a flow diagram of the method of the present invention. Thus, the following method 20 is implemented in the base unit 3 of the OPS imaging system. The method 20 is based on a physical model that relates the intensity distribution within the focal plane F to the quantity of light scattered in the region between the focal plane F (where vessels reside) and the tissue surface 12A.

The first step 22, requires that the polarized light 11 be emitted at the target tissue 12. In step 24, after the AOP 10 removes the surface-reflected light 13A, the de-polarized light 13B is passed to the base unit 3 where the measured light intensity $I^{M}(x)$ forms an "original" or "subsampled" image. In step 26, the processor 4 estimates the intensity of the light scattered in the foreground tissue for every pixel 36, $I_{s}(x)$. In step 28, the processor 4 subtracts the scattered light $I_{s}(x)$ from the original image $I^{M}(x)$ to determine light intensity in the focal plane F, $I_{r}(x)$. In step 30, using this "corrected light intensity," Beer's Law is then used to determine the concentration of hemoglobin.

In particular, in estimating the light scattered (I_s) in the foreground tissue 32 for every pixel 36 position, step 26 comprises:

a) Calculating the probability (P_s) that a photon originating from position x_i (x_i ≠
 x) on the focal plane F scatters within foreground tissue 32 into the optical

path 38 of the pixel 36 at x on the tissue surface 12. The probability (P_s) can be estimated using standard Monte Carlo computation methods.

b) Centering the probability function (P_s) on x and multiplying P_s with the image intensity distribution I^M . Summing over all x_i in the focal plane F to obtain the total intensity of scattered light I_s at x:

$$I_s(x) = \int_{image} (dx_i)^2 P_s(|x - x_i|, z_f) \cdot I^M(x_i)$$

Step 28, subtracting scattered light $I_s(x)$ from the original image $I^M(x)$ to determine intensity in focal plane F, $I_t(x)$, is given by:

10
$$I_f(x) = \frac{I^M(x) - I_s(x)}{\exp(-\mu_t z_f)}.$$

5

Next, step 30, using Beer's Law to determine the concentration of hemoglobin [Hb] is given by:

$$[Hb] = -\frac{1}{\alpha D} \log \left(\frac{I_f(x_v)}{I_f(x_b)} \right) = -\frac{1}{\alpha D} \log \left(\frac{I^M(x_v) - I_s(x_v)}{I^M(x_b) - I_s(x_b)} \right),$$

where x_{ν} and x_{b} are position of vessel and background, respectively; D is the diameter of

the vessel BV and α is the absorptivity of Hb.

5 .

10

15

The parameters which contribute to the light scattered from other positions on the focal plane $F(I_s(x))$ are (1) depth z_r (parameter of probability function); (2) scattering coefficient (also a parameter of the probability function); and (3) light distribution in the focal plane F(e.g., shadows, diameter of main vessel).

The method 20 of the present invention can be implemented in a computer program for use by the processor 4. As shown in the Appendix, by way of example only, computer code is written in MatLab language to perform the method 20 of the present invention. Fig. 4 provides an overview of the flow of the computer program. In particular, step 24 comprises generating the subsampled or original image from the measured light intensity, $I^{M}(x)$. Step 26 comprises sub-steps 26A-26C, namely, calculating the scattering probability function (P_s) 26A, convolving the subsampled image and the probability function P_s 26B and generating the scattered light image, $I_s(x)$ 26C. Next, step 28 comprises sub-steps 28A-28B: in step 28A, a scattered light image is subtracted from the subsampled image (forming step 28A) and a scattering corrected image based on $I_r(x)$ is generated in step 28B. Although not shown in the computer program attached as Appendix, the processor 4 then calculates the concentration of hemoglobin [Hb] using Beer's Law, discussed earlier.

APPENDIX

```
t0 = clock;
              Lpixel = 0.0002; %cm
 5
              load ImgScat.mat;
              fprintf ('Assumed focal depth = %.0f
       microns/n', z_focus*10^4);
              mu t = mu a + mu sp;
              %% Resize image so that pixel size is same as for ImgScat.
              ImgSize = size(Img);
10
              subsamp = dx im/Lpixel;
              NewSize = ImgSize/subsamp;
              if ~ (fix(NewSize) ==NewSize)
                 fprintf ('Sampling interval incorrect. \n')
15
                 return;
              end
              Img2 = imresize (Img, NewSize, 'bilinear');
              %% Determine intensity at focal plane.
              Focalimg = Img2/(exp(-mu t*z focus) + TotCylProb);
              if ('simulate' == simulate)
20
                Focallmg = Img2;
              end
              ImgScat = [ImgScat (:, end:-1:2) ImgScat];
              ImgScat = [ImgScat (end:-1:2,:); ImgScat];
              %% Determine amount of scattered light in foreground.
25
              Ifore = conv2 (FocalImg, ImgScat, 'same');
              Ifore = imresize (Ifore, ImgSize, 'bilinear');
              ImCorr = Img - Ifore;
              if ('simulate' == simulate)
30
                 ImgCorr = Img*exp(-mu t*z focus) + Ifore;
              end
              fprintf (' \lg \max = \%.3f \n', \max (Ifore(:)));
              fprintf (' lg min = \%.3f \n', min (Ifore(:)));
```

function [ImgCorr, Ifore] = ForeScat (Img, simulate);

%%% Calculate the amount of light scattered within the foreground tissue.

```
fprintf (' lg mean = %.3f \n', mean (lfore (:)));
              fprintf ('Time2 = \%.1f \n', etime (clock, t0) );
       return;
       function [rho, prob, CylProb] = ScatProbAll ( z_focus, mu_a,
 5
       mu sp, twoD);
       %%% Calculate the probabilities of foreground scattering from
       all pixels in
       %%% focal plane. This function calls ScatProbl for all pixels.
       %%% Robert Danen (5.21.99)
10
              t0 = clock:
       %%
              Definitions ----
              Lpixel = 0.0002; %cm = size of pixel
              NA = 0.165; %numerical aperture of probe
              NA theta = asin (NA);
15
              TanNA = tan ( NA theta );
              sNA = 2*pi* (1-cos(NA theta)); %solid angle of NA
              Setup sampling intervals -----
       %%
              if (twoD)
                    ScatRegion = 0.024; %0.07; %cm
                    subsamp = 2; %4;
20
              dx im = subsamp*Lpixel; %sampling of final probability
       matrix
                    dx = dx im/4;
                    dx 2 = dx/2;
25
              else
                    ScatRegion = 0.07; %0.128 %cm
                    dx im = 0.0006; %cm
                    dx = dx im;
30
              rho = [0:dx:ScatRegion]; %rho = [0:dx:0.1278];
       %%
              Do calculations -----
              mu_t = mu_a + mu_sp;
              clear prob;
              prob(1) = 0;
              for i=2:length(rho) %DO NOT calculate for rho=0.
35
                    rho1 = rho(i);
                    prob(i) = ScatProb1 (rho1, z focus, TanNA, Lpixel, mu_a,
       mu_sp);
```

```
%prob = ScatProbV(rho1, 0.0025, z focus, 1.25*z focus,
       250, 57); prob(i) = sum(prob);
              end
              DetNormA = ( (dx/Lpixel)^2) /sNA; %Normalize intensity
              DetNormB = ( (dx_im/Lpixed)^2) /sNA; %Normalize intensity
 5
       %%-
              IN-VITRO FACTOR
       %%
              prob = 5.0*prob;
              prob = 2.7*prob; %for g=0.8
       %
       %%-
10
              Uniform (cylindrical) geometry %%
       %%
              CylProb = DetNormA*(2*pi/dx)*rho.*prob;
              TotCylProb = sum(CylProb);
              fprintf ('Scattering factor = %.4f \n', TotCylProb);
15
       %% Make 2D matrix containing scattering probabilities; save to
       file.
              If(twoD)
                     [X, Y] = meshgrid (0:dx im:rho(end), 0:dx im:rho(end));
                     Z = sqrt(X.^2+Y.^2);
20
                     [sv, sx] = size(Z);
                     fprintf (' Quarter-size of probability image - %d x
       %d\n',sy,sx);
              ImgScat = zeros(sy,sx);
              for iy=1:sy
25
                     rhol = Z(iy,ix);
                     ind = find ( rho >=(rho1-dx 2) & rho < (rho1+dx 2) );
                     if (~isempty(ind))
                           ImgScat (iy,ix) = prob(ind) *DetNormB;
                     else
30
                           ImgScat(iy,ix) = 0.0;
                     end
              end
       end
       save 'c:\RadiativeTransfer\ImgScat' dx im z focus mu a
35
       mu sp TotCylProb ImgScat;
              %save 'c:\RadiativeTransfer\ImgScat.m' ImgScat -ASCII;
       end
       fprintf('Time1 = \%.1f \n', etime(clock,t0));
```

```
function prob = ScatProb1(rho1, z_focus, TanNA, Lpixel, mu_a,
       mu sp)
       %% Calculates probability that a photon originating at a
 5
       %% rho1 in the focal plane from the pixel of interest scatters
       into
       %% the collection cone of the pixel.
       %% Robert Danen (5.21.99)
              t0 = clock;
10
              g = 0.8;
       %%
              Definitions----
              mu_t = mu_a + mu_sp;
              Dtheta2 = 0.01; %radians
              Dphi3 = 0.01; %radians
              Drho3 = 0.00002; %cm
15
       %%
              Do calculations----
              %Theta2_lim = atan(z_focus/rho1);
              Theta2 \lim = abs (atan(1.0/((rho1/z focus) - TanNA)));
              Theta2 = [Dtheta2:Dtheta2:Theta2_lim];
20
              CosTheta2 = cos(Theta2); %cos theta
              SinTheta2 = sin(Theta2);
              r2 = rho1./CosTheta2:
              SA source = 2*sin(0.5*Dtheta2)*Dphi3;
              source = SinTheta2.*SA_source;
25
              clear Zprob;
              for ith2 = 1:length(Theta2)
                     %fprintf('ith2= %d th2=
       %.3f\n'ith2,acos(CosTheta2(ith2));
              phi3_lim_a = atan( SinTheta2(ith2) * TanNA );
              phi3 lim = fix(phi3 lim a/Dphi3)*Dphi3;
30
              if (phi3 \lim == phi3 \lim a)
                phi3 1im = phi3 lim - 1;
              end
              phi3 = pi + [0:Dphi3:phi3 lim];
35
              %% determine length where scattering event can occur
       (from rho3 lim1 to rho3 lim2)
```

```
[rho3 lim1, rho3 lim2] = ConeIntersection(phi3,
        CosTheta2 (ith2), SinTheta2 (ith2), rho1, TanNA);
               r1 lim2 = CoordinateTrans (rho3 lim2, phi3,
        CosTheta2 (ith2), SinTheta2(ith2), rho1);
               ind = find(r1_lim2(:,3) > z_focus | r1_lim2(:,3) < 0);
 5
               rho3 lim2 (ind) = -z focus./(
        SinTheta2(ith2)*cos(phi3 (ind));
               %%% Integrate over rho3 for each phi3. %%%
               clear Zprob1:
               for iphi3 = [1:length(phi3)]
10
                     clear rho3;
                     rho3 = [rho3 \ lim1 \ (iphi3):Drho3:rho3 \ lim2 \ (iphi3)];
                      %clear p1 p2 p3;
                     p1 = exp(-mu_t^*rho3);
15
                     r1 s = CoordinateTrans(rho3, phi3(iphi3),
        CosTheta2(ith2), SinTheta2 (ith2), rho1);
               %r1 \text{ out} = z \text{ focus } -r1 \text{ s(:,3)};
               r1 s M2 = sum(r1 s.^2,2);
               r1 s M2 = r1 s M2';
               r1 s M = sqrt(r1 s M2);
20
               SecTheta1 = r1 s M./(r1 s(:,3));
               r1_out = (z_focus - r1_s(:,3) ').*SecTheta1;
               p3 = exp(-mu t*r1 out);
               dtheta1 = 0.5*Lpixel./(r1 s(:,3)');
25
               ind = find(dtheta1 > TanNA);
               dtheta1(ind) = TanNA;
               dtheta1 = atan( dtheta1 );
               SApixel = 4*dtheta1.*sin(dtheta1); %=
        4*dtheta1.*sin(dtheta1)
               %theta s = acos((rho3.^2 + r1 s M2 - r1))
30
        rho1^2)./(2*r1 s_m.*rho3));
               %ptheta = SApixel.*HGPhase(theta s,g);
               ptheta = Sapixel/(4*pi); % (for mu s') % =)(solid
        angle of pixel)/4pi
35
               Zprob1(iphi3) = sum(ptheta.*p1.*p3);
               Zprob(ith2) = 2*sum(Zprob1) - Zprob1(1); %Integrate over
        only half the NA cone to save computational time.
```

```
end
```

```
Zprob = source.*Zprob;
            p2 = 1 - exp(-mu sp*Drho3);
            prob = p2 * sum( Zprob ); %multiply by 2 because only
 5
      integrated over half of cone
       %
            fprintf('elapsed time = %.1f \n', etime(clock,t0) );
       %
            fprintf('prob = \%.5e \n', prob);
      return;
       10
      function [rho3 a, rho3 b] = ConeIntersection(phi3, CosTheta2,
       SinTheta2, xshift1, TanNA)
            Calculate roots of cone in cylindrical frame
       %%
            phi3 is a vector, remaining inputs must be scalar.
       %%
            CosPhi3
                        = \cos (phi3);
            SinPhi3
15
                        = \sin (phi3);
            CosPhi3 2
                        = CosPhi3.^2;
            SinPhi3 2
                        = SinPhi3.^2;
            CosTheta2 2 = CosTheta2^2;
            SinTheta2_2 = SinTheta2^2;
20
            TanNA .2
                         = TanNA<sup>2</sup>;
                   = (SinTheta2 2 * TanNA 2) * CosPhi3 2);
            denom = (term1 - CosTheta2 2*CosPhi3 2 - SinPhi3 2);
            term2 = sqrt (term1 - SinPhi3 2);
25
            term3 = CosTheta2*CosPhi3;
            term4 = xshift1./denom;
            rho3 a = \text{term4.*(term3 + term2)};
            rho3 b = term4.*(term3 - term2);
      return;
30
      function r1 = CoordinateTrans (rho3, phi3, CosTheta2,
       SinTheta2, xshift1)
       %% Transform from source to detection coordinate frames, for
      z3=0.
       %% rho and phi3 are vectors with same index, remaining inputs
35
```

```
must be scalar.
```

```
CosPhi3 = cos(phi3);

SinPhi3 = sin(phi3);

r1(1,:) = CosTheta2*(rho3.*CosPhi3) + xshift1;

r1(2,:) = rho3.*SinPhi3;

r1(3,:) = -SinTheta2*(rho3.*CosPhi3);

r1 = r1':
```

return;

function prob = HGPase (theta,g);
%%% Normalized Henyey-Greenstein phase function.

$$a = 1 - g^2;$$

 $b = 2*(1 + g^2 - 2*g*cos(theta)).^1.5;$
 $prob = a./(2*pi*b);$

15 return;

10